



APRONS

Revolution Front Elastic Back Base

B403



Shown in SH-18 Shape Camo

- Elastic back panels for increased comfort & support
- Classic Nylon Colors
- Professionally designed for an exceptional fit and comfort
- SmartTrack enabled

Model Number	Description
B403	Revolution Front Elastic Back Base
Colors	Classic Nylon Colors
Core	Lightweight Lead KIARMOR Lead-Free
Protection	.25mm LE, .35 mm LE, .50mm LE
Sizes	Male S-5XL Female XS,-5XL Length alterations $\pm 2"$ or $\pm 4"$
Fabric	Nylon
Monogram	Block or <i>Script</i>



REV002

1040 Avenida Acaso
Camarillo, CA 93012, USA

tel: 805.987.5255
fax: 805.482.8424

info@infabcorp.com
www.infabcorp.com

