



APRONS

Revolution Full Wrap Black Belt

203




Shown in 506 Teal Black - Block

- Partial overlap to increase protection
- Hook and loop closure side slits for increased mobility
- Cool Wear moisture-wicking interior fabric
- Professionally designed for an exceptional fit and comfort
- Snap on, fully adjustable Black Belt®
- SmartTrack enabled

Model Number	Description
203	Revolution Full Wrap Black Belt
Colors	500, 600, 700, 800, 900, Remix
Core	Lightweight Lead KIARMOR Lead-Free
Protection	.25mm LE, .35mm LE, .50mm LE
Sizes	Male S-5XL Female XS-5XL Length alterations $\pm 2"$ or $\pm 4"$
Fabric	Revolution Fabric with Cool Wear
Optional Base	B203 Classic Nylon Colors without Cool Wear (Custom alterations)
Monogram	Block or <i>Script</i>

Shown in 602 Black/Fuschia

Hook and loop
closure side slits

REV002

1040 Avenida Acaso
Camarillo, CA 93012, USAtel: 805.987.5255
fax: 805.482.8424info@infabcorp.com
www.infabcorp.com