



BARRIERS

Wide Mobile Barrier

683482



- Powder coated steel frames, precision TIG welded for years of trouble free service completely enclosing the leaded glass panel for increased durability and safety
- Lead glass
- 1.5 mm lead equivalent opaque base covered in "Sintra" plastic for easy cleaning and durability
- Modular parts can be easily replaced, including the glass if needed.

Model Number	Description
683482	Wide Mobile Barrier - Lead
Dimensions	72" H x 52½" W X 25" D
Weight	235 lbs.
Window	36" H x 48" W
Glass	2.0 mm Pb LE
Frame Color	White
Frame	Powder Coated Steel
Curtain Color	Revolution Grey
Wheels	4 swivel casters, 2 with brakes
Crate Size	30" x 72" x 60"
Crate Weight	480 lbs / 217.72 kg



REV002

1040 Avenida Acaso
Camarillo, CA 93012, USA

tel: 805.987.5255
fax: 805.482.8424

info@infabcorp.com
www.infabcorp.com

