



Dental X-Ray Protection

Infab Corporation

3651 Via Pescador
Camarillo, CA, USA 93012
Phone: 805-987-5255
Fax: 805-482-8424

Contact Us

M-F 6:30am-5pm PDT
info@infabcorp.com
www.InfabCorp.com



Custom Apron Options

Custom Logos

Use your logo on your apron.



Monograms

Block
Linda Joyce, MD

Script
Linda Joyce, MD

Colors

Outside, Inside, and Trim Options



Over 2,000 Colors at
Infabcorp.com/colors



COLORS

Protection



SMART ARMOR

10 YEAR CORE GUARANTEE



Regular Lead

Composed of emulsion polymers, finely divided 100% pure lead metal particles, stabilizers and pigments.



Lightweight Composite

Composed of emulsion polymers, finely divided metal particles including lead, antimony and barium, stabilizers and pigments.



Greenlite Lead-Free

Composed of emulsion polymers, finely divided lead-free metal particles including antimony and barium, stabilizers and pigments.



REVOLUTION

Introducing a Revolution in Protection

Lead & Toxin-Free, Super Flexible, Superior Radiation Protection at 60 - 150 kVp
20% Lighter than Greenlite. Up to 50% Lighter than any competitor's product.

10 Year Warranty

Infab guarantees our core material will not stiffen, crack, or crumble for a minimum of 10 years from the date of purchase.

Latest Technology

With over 30 years of research, Infab's 3 types of Smart Armor core material are the best in the industry.

Highest Attenuation

Infab's Smart Armor line (Greenlite lead-free, Lightweight Composite, and Regular Lead) of core material have the highest attenuation level in the industry from 60 kVp to 130 kVp direct beam

Adult

Adult Dental Apron



FEATURES

- 24" W x 27" L
- Non-slip inner material
- .3mm Pb equivalency

MODEL NUMBER

MODEL	ONE SIZE FITS ALL (Custom Sizes Available)
DTL-ADULT	24" W x 27" L

Adult Dental with Collar



- 24" W x 27" L
- Genuine VELCRO® brand closure around neck
- Non-slip inner material
- .3mm Pb equivalency
- Attached thyroid protection

MODEL NUMBER

MODEL	ONE SIZE FITS ALL (Custom Sizes Available)
DTL-ADULT-TC	24" W x 27" L

Pano Cape Apron



- 23" W x 25" L
- Split in back with Genuine VELCRO® brand closure
- .3 mm Pb equivalency

MODEL NUMBER

MODEL	ONE SIZE FITS ALL (Custom Sizes Available)
DTL-PANO	23" W x 25" L

Choose a Lead Type:



Choose Color Options



Over 2,000 Colors at
Infabcorp.com/colors

Technician



Coat Style Technician

- 24" W x 36" L
- Arm holes
- Tie straps
- .3mm Pb equivalency

MODEL NUMBER

MODEL	ONE SIZE FITS ALL <i>(Custom Sizes Available)</i>
--------------	---

DTL-TECH	24" W x 36" L
-----------------	---------------

Childrens

Children's Dental Apron



- 18" W x 23" L
- Non-slip inner material
- .3mm Pb equivalency

MODEL NUMBER

MODEL	ONE SIZE FITS ALL <i>(Custom Sizes Available)</i>
--------------	---

DTL-CHILD	18" W x 23" L
------------------	---------------

Children's Dental with Collar



- 18" W x 24" L
- Genuine VELCRO® brand closure around neck
- Non-slip inner material
- .3mm Pb equivalency
- Attached thyroid protection

MODEL NUMBER

MODEL	ONE SIZE FITS ALL <i>(Custom Sizes Available)</i>
--------------	---

DTL-CHILD-TC	18" W x 24" L
---------------------	---------------

Eye & Thyroid



Single Use Eye & Thyroid Shields


Infab Corporation is pleased to introduce patient protection eye and thyroid shields to our line of revolutionary Smart Armor protection products. Thoughtfully and carefully engineered to maximize patient protection and comfort with no noise or image artifacts. Highly effective during CT scans, dental radiography or any procedure where X-Ray radiation is present.

Increase patient comfort and confidence

Recent news stories focusing on the harmful effects to the eyes and thyroid from CT and dental scans can make patients uneasy. Now there's a simple, safe and effective way to offer shielding without compromising image quality.



Eye Shields



Back




Minimum 40% Dose Reduction

- Non Latex
- No heavy metals
- Sanitary single use
- Disposable
- No Noise, specially designed to avoid artifacts in the image
- Package of 50
- Adhesive backing to hold in place and for ease of positioning

Model	Description
PP-ES	Smart Armor Eye Shield - Box of 50

Thyroid Shield



Minimum 60% Dose Reduction

- Non Latex
- No heavy metals
- Sanitary single use
- Disposable
- No Noise, specially designed to avoid artifacts in the image
- Package of 50

Model	Description
PP-TS	Smart Armor Thyroid Shield - Box of 50



Barium sulfate shields are non-latex and 100% disposable

