

Classic Thinking Cap



- Head and Skull Radiation Protection
- Available in KIARMOR Bi-layer Lead-Free, KIARMOR Bi-layer Low-Lead, Lead
- Classic nylon color options
- .25 mm LE, .35 mm LE, .5 mm LE
- Inside lined with Cool Wear

SIZE: Small, Medium, Large, X-Large, XX-Large
COLORS: Classic Nylon Colors
PROTECTION: .25 mm LE, .35 mm LE, .50 mm LE
CORE: Lead, KIARMOR Bi-Layer Low-Lead, KIARMOR Bi-Layer Lead-Free

Size	Description	PROTECTION OPTIONS		
				
Small	.25 mm LE Thinking Cap - 21" (53.34 cm)	399998-RL	399998-KI-LF	399998-KI-LL
Medium	.25 mm LE Thinking Cap - 22" (55.88 cm)	399999-RL	399999-KI-LF	399999-KI-LL
Large	.25 mm LE Thinking Cap - 23" (58.42 cm)	400000-RL	400000-KI-LF	400002-KI-LL
X-Large	.25 mm LE Thinking Cap - 24" (60.96 cm)	399991-RL	399991-KI-LF	399991-KI-LL
XX-Large	.25 mm LE Thinking Cap - 25" (63.50 cm)	399992-RL	399992-KI-LF	399992-KI-LL
Medium	.50 mm LE Thinking Cap - 22" (55.88 cm)	400002-RL	400002-KI-LF	400002-KI-LL
Large	.50 mm LE Thinking Cap - 23" (58.42 cm)	400003-RL	400003-KI-LF	400003-KI-LL
X-Large	.50 mm LE Thinking Cap - 24" (60.96 cm)	399993-RL	399993-KI-LF	399993-KI-LL
XX-Large	.50 mm LE Thinking Cap - 25" (63.50 cm)	399994-RL	399994-KI-LF	399994-KI-LL