



ACCESSORIES

Revolution Hook and Loop Closure Back Thinking Cap

400025



Model	Description
400025	Revolution Hook and Loop Closure Back Thinking Cap
Features	Exclusive Revolution outer fabric Scatter radiation protection Cool Wear Moisture-Wicking Fabric Hook & loop back closure for customizable fit
Sizes	One Size Fits All
Colors	700 Series Revolution Colors
Core	Lightweight Lead KIARMOR Bi-Layer Lead-Free
Protection	.25mm LE, .35 mm LE, .50mm LE
Fabric	Revolution Fabric
Monogram	Block or <i>Script</i>