

# Shin Guards



Model	Description
444000	Small Shin Guard 11" W x 12" H (Not including foot protection)
444005	Medium Shin Guard 12" W x 13" H (Not including foot protection)
444010	Large Shin Guard 14" W x 14" H (Not including foot protection)

**COLORS:** Classic Nylon Colors  
**SIZES:** Small, Medium, Large  
**CORE:** Lightweight Lead, KIARMOR Bi-Layer Low-Lead, KIARMOR Bi-Layer Lead-Free

- Reduce radiation exposure to the shins
- Hook and loop closures for easy adjustments