

Optional Right or Left Sleeve SLV



Colors	Revolution 700 Series & All Classic Nylon Colors
Attached	Sewn on, snaps or hook and loop closure
Core	Lightweight Lead, KIARMOR Bi-Layer Low-Lead, KIARMOR Bi-Layer Lead-Free

Model	Description	
SLV-L	Left Sleeve	Sewn On
SLV-L-S	Left Sleeve	Attached by Snaps
SLV-L-V	Left Sleeve	Attached by hook and loop closure
SLV-R	Right Sleeve	Sewn On
SLV-R-S	Right Sleeve	Attached by Snaps
SLV-R-V	Right Sleeve	Attached by hook and loop closure

- Added protection of the arm and shoulder
- .50mm LE
- Available on most Infab apron models
- Optional removable snap on, permanent sewn on or hook and loop closure attached configuration