


**10 Arm Budget Saver Mobile Apron Rack
683427**



- Economically priced
- Heavy gauge steel
- Powder coated
- 2 high-quality hospital grade locking casters
- Stores up to 10 protection aprons
- Convenient additional storage on base
- Optional vest and skirt adapter (Part Number 683000)

Model Number	Description
683427	Budget Saver Mobile Apron Rack with 10 arms
Dimensions	58" H x 44½" W x 30" D
Colors	White, Custom Options
Options	Powder Coated Steel Stainless Steel (Upgrade)
Wheels	4 swivel casters, 2 with brakes
Shipping	Box 1: 51" x 7" x 7", 36 lbs / 16.33 kg Box 2: 16" x 32" x 13", 34 lbs / 15.42 kg Box 3: 58" x 16" x 9", 40 lbs / 18.14 kg Box 4: 22" x 44" x 2", 21 lbs / 9.53 k