


20 Hanger Mobile Apron Valet

683438



- Space saving; providing major storage while taking up minimum space.
- Steel construction with baked enamel finish.
- 2 high-quality hospital grade locking casters
- 20 Deluxe Custom Apron Hangers included
- Optional: full circle loop hangers (-H)

Model Number	Description
683438	Mobile Apron Valet with 20 Hangers
683450	Additional hangers (ea)
683451	Additional hangers (set of 5)
<hr/>	
Dimensions	58" H x 48½" W X 25" D
Colors	White, Custom Options
Options	Powder Coated Steel Stainless Steel (Upgrade)
Wheels	4 swivel casters, 2 with brakes
Shipping	Box 1: 55" x 7" x 7", 32 lbs / 14.51 kg Box 2: 16" x 32" x 14", 43 lbs / 19.50 kg