


**5 Arm Budget Saver Mobile Apron Rack
683426**



- Economically priced
- Heavy gauge steel
- Powder coated
- 2 high-quality hospital grade locking casters
- Stores up to 5 protection aprons
- Convenient additional storage on base
- Optional vest and skirt adapter (Part Number 683000)

Model Number	Description
683426	Budget Saver Mobile Apron Rack with 5 arms
Dimensions	58" H x 26" W X 30" D
Colors	White, Custom Options
Options	Powder Coated Steel Stainless Steel (Upgrade)
Wheels	4 swivel casters, 2 with brakes
Shipping	Box 1: 51" x 7" x 7", 31 lbs / 14.06 kg Box 2: 16" x 32" x 14", 34 lbs / 15.42 kg Box 3: 16" x 41" x 7", 21lbs / 9.53 kg Box 4: 26" x 21" x 2", 12 lbs / 5.44 kg