



# Wall Mounted Apron Rack

## 683429



- Heavy gauge steel
- 200 lbs. Capacity
- Powder coated
- Scratch resistant
- Easy to install
- Requires only minimal space
- Stores up to 10 aprons
- Optional vest and skirt adapter (Part Number 683000)

| Model Number      | Description                                      |
|-------------------|--|
| <b>683429</b>     | Wall Mounted Apron Rack                          |
| <b>Dimensions</b> | 12" H x 36" W X 26" D                            |
| <b>Colors</b>     | White, Custom Options                            |
| <b>Options</b>    | Powder Coated Steel<br>Stainless Steel (Upgrade) |
| <b>Shipping</b>   | Box 1: 40 lbs (16" x 58" x 9")                   |