


# 2" Wide Belt WBT



- Comfortable shoulder pads
- Criss-cross straps with 2" belt buckles
- Double welt pockets



Model Number	Description
WBT	2" Wide Belt
Colors	Classic Nylon Colors
Core	Lightweight Lead KIARMOR Bi-Layer Low Lead KIARMOR Bi-Layer Lead-Free
Protection	.25mm LE, .35 mm LE, .50mm LE
Sizes	Male S-5XL Female XS-5XL
Monogram	<b>Block</b> or <i>Script</i>

Model	Description	Gender	Sizes	Color	Protection Options		
					Lightweight Lead	KIARMOR Low Lead	KIARMOR Lead-Free
WBT	2" Wide Belt	M	S, M, L, XL, 2XL, 3XL, 4XL, 5XL	Classic Nylon Colors	LL	KI-LL	KI-LF
		F	XS, S, M, L, XL, 2XL, 3XL, 4XL, 5XL	Classic Nylon Colors	LL	KI-LL	KI-LF