



IN STOCK & READY TO SHIP!

Model Number:

UFP • Lightning Front Protection

Key Features:

Full front panel for maximum protection that distributes weight between waist and shoulder with hook and loop closures and a double welt pocket.

Colors:





Charcoal (205) Blue (108) Red (105) Purple (119) Teal (117) Black (114)

Unisex Sizing:

Size	Chest	Waist	Length
S	34-37" / 86-94 cm	28-31" / 71-79 cm	38" / 97 cm
M	38-41" / 97-104 cm	32-35" / 81-89 cm	38" / 97 cm
L	42-45" / 107-114 cm	36-39" / 91-99 cm	38" / 97 cm
XL	46-49" / 117-125 cm	40-43" / 102-109 cm	42" / 107 cm
2XL	50-52" / 127-132 cm	44-46" / 112-117 cm	44" / 112 cm

Lead Options:

CORE	STANDARD RADIATION PROTECTION CONFIGURATION		
	PROTECTION	.35MM LE	.50MM LE
	ATTENUATION	98.0%	99.4%
	ABSORBED DOSE	2.0%	0.6%
	PROTECTION	.35MM LE	.50MM LE
	ATTENUATION	97.6%	99.1%
	ABSORBED DOSE	2.45%	0.9%

Add On:

UST-S1 • Thyroid Collar without Binding

Key Features:

A comfortable fit that is fully adjustable with hook and loop closure

* Color options are the same listed above.

