



## IN STOCK & READY TO SHIP!

Model Number:

## UVS • Lightning Vest & Skirt

### Key Features:

Partial overlapping front panels for maximum protection with hook and loop closures, double welt pocket, and a belt with buckle closure on skirts.

### Colors:



Charcoal (205) Blue (108) Red (105) Purple (119) Teal (117) Black (114)

### Unisex Sizing:

Size	Chest	Waist	Vest Front	Skirt
S	34-37" / 86-94 cm	28-31" / 71-79 cm	24" / 61 cm	20" / 51 cm
M	38-41" / 97-104 cm	32-35" / 81-89 cm	24" / 61 cm	22" / 56 cm
L	42-45" / 107-114 cm	36-39" / 91-99 cm	24" / 61 cm	23" / 58 cm
XL	46-49" / 117-125 cm	40-43" / 102-109 cm	24" / 61 cm	24" / 61 cm
2XL	50-52" / 127-132 cm	44-46" / 112-117 cm	26" / 66 cm	25" / 64 cm

### Lead Options:

CORE	STANDARD RADIATION PROTECTION CONFIGURATION		
	PROTECTION	.35MM LE	.50MM LE
	ATTENUATION	98.0%	99.4%
	ABSORBED DOSE	2.0%	0.6%
	PROTECTION	.35MM LE	.50MM LE
	ATTENUATION	97.6%	99.1%
	ABSORBED DOSE	2.45%	0.9%

\*BACK PROTECTION ON BOTH OPTIONS: 0.25MM LE

### Add On:

## UST-S1 • Thyroid Collar without Binding

### Key Features:

A comfortable fit that is fully adjustable with hook and loop closure

\* Color options are the same listed above.

