

Revolution Full Wrap Elastic Back

603

- Cool Wear moisture-wicking interior fabric
- Professionally designed for an exceptional fit and comfort
- SmartTrack enabled



Shown in 616 Black/Volt

Model Number	Description
603	Revolution Full Wrap Elastic Back
Colors	500, 600, 700, 800, 900, Remix
Core	Lightweight Lead KIARMOR Bi-Layer Low Lead KIARMOR Bi-Layer Lead-Free
Protection	.25mm LE, .35 mm LE, .50mm LE
Sizes	Unisex XS-2XL Length Alterations $\pm 2''$ or $\pm 4''$
Fabric	Revolution Fabric with Cool Wear
Optional Base	B603 Classic Nylon Colors without Cool Wear (Custom alterations)
Monogram	Block or <i>Script</i>



Sizing

	Bust	Waist	Length
X-Small	30-33" / 76-84 cm	25-27" / 64-69 cm	38" / 97 cm
Small	34-37" / 86-94 cm	28-31" / 71-79 cm	38" / 97 cm
Medium	38-41" / 97-104 cm	32-35" / 81-89 cm	38" / 97 cm
Large	42-45" / 107-114 cm	36-39" / 91-99 cm	38" / 97 cm
X-Large	46-49" / 117-125 cm	40-43" / 102-109 cm	38" / 97 cm
XX-Large	50-52" / 127-132 cm	44-46" / 112-117 cm	38" / 97 cm

Radiation Protection Options

STANDARD CONFIGURATION		OPTION 1		OPTION 2	
.50mm LE		.35mm LE		.25mm LE	
Absorbed Dose	Attenuation	Absorbed Dose	Attenuation	Absorbed Dose	Attenuation
0.9%	99.1%	2.45%	97.6%	4.9%	95.1%
0.6%	99.4%	1.9%	98.1%	3.5%	96.5%
0.6%	99.4%	2.0%	98.0%	3.9%	96.1%
Highest Protection		Best Combination of Protection & Weight		Lightest Weight	



Model	Description	Gender	Sizes	Color	Protection Options		
					Lightweight Lead	KIARMOR Low Lead	KIARMOR Lead-Free
603	Revolution Full Wrap Elastic Back	UNISEX	S, M, L, XL, 2XL, 3XL, 4XL, 5XL	Revolution	LL	KI-LL	KI-LF