

Revolution Full Overlap Lumbar Vest & Skirt

703

The 703 brings the Classic Infab design into the Revolution line, adding a Black Belt for additional lumbar support. This full overlap apron is lightweight, offers great protection, and has a Classic fit.

- Full overlap
- Cool Wear moisture-wicking interior fabric
- Professionally designed for an exceptional fit and comfort
- Built-in Black Belt® in vest and skirt for better lumbar support
- SmartTrack enabled



Shown in 505 Navy/Black



KTC Buckle to keep belt attached and in place.



All new color block design for 703. Inside panel of vest & skirt

Model Number	Description
703	Revolution Full Overlap Lumbar Vest & Skirt
Colors	500, 600, 700, 800, 900, Remix
Core	Lightweight Lead KIARMOR Bi-Layer Low Lead KIARMOR Bi-Layer Lead-Free
Protection	.50mm LE
Sizes	Male S-5XL Female XS-5XL Length alterations ±2" or ±4"
Fabric	Revolution Fabric with Cool Wear
Monogram	Block or <i>Script</i>

STANDARD RADIATION PROTECTION CONFIGURATION	
.50mm LE	
Absorbed Dose	Attenuation
0.9%	99.1%
0.6%	99.4%
0.6%	99.4%
Highest Protection	



Model	Description	Gender	Sizes	Color	Protection Options		
					Lightweight Lead	KIARMOR Low Lead	KIARMOR Lead-Free
703	Revolution Full Overlap Lumbar Vest & Skirt	M	S, M, L, XL, 2XL, 3XL, 4XL, 5XL	Revolution	LL	KI-LL	KI-LF
		F	XS, S, M, L, XL, 2XL, 3XL, 4XL, 5XL	Revolution	LL	KI-LL	KI-LF