


Revolution Lumbar Vest & Skirt Base

BL103



- “Back Flap” for better weight balance
- Built-in “Black Belt” in Vest and Skirt for better lumbar support
- Partial overlap to increase protection
- Professionally designed for an exceptional fit and comfort
- Smart Track enabled

Model Number	Description
BL103	Revolution Lumbar Vest & Skirt Base
Colors	Classic Nylon Colors
Core	Lightweight Lead KIARMOR Bi-Layer Low Lead KIARMOR Bi-Layer Lead-Free
Protection	.25mm LE, .35 mm LE, .50mm LE
Sizes	Male S-5XL Female XS-5XL Length alterations ±2”or ±4”
Fabric	Nylon
Monogram	Block or <i>Script</i>

Radiation Protection Options

STANDARD CONFIGURATION		OPTION 1		OPTION 2	
.50mm LE		.35mm LE		.25mm LE	
Absorbed Dose	Attenuation	Absorbed Dose	Attenuation	Absorbed Dose	Attenuation
0.9%	99.1%	2.45%	97.6%	4.9%	95.1%
0.6%	99.4%	1.9%	98.1%	3.5%	96.5%
0.6%	99.4%	2.0%	98.0%	3.9%	96.1%
Highest Protection		Best Combination of Protection & Weight		Lightest Weight	



Back panel is always .25mm LE unless requested otherwise

Shown in SH-23 Art Deco

Model	Description	Gender	Sizes	Color	Protection Options		
					Lightweight Lead	KIARMOR Low Lead	KIARMOR Lead-Free
BL103	Revolution Lumbar Vest & Skirt Base	M	S, M, L, XL, 2XL, 3XL, 4XL, 5XL	Classic Nylon Colors	LL	KI-LL	KI-LF
		F	XS, S, M, L, XL, 2XL, 3XL, 4XL, 5XL	Classic Nylon Colors	LL	KI-LL	KI-LF