


Economy Mobile Curtain

REV-MC



A unique new design for operating physicians that allows unparalleled protection for close-to-table studies. Sturdy, yet thoughtfully engineered with a small footprint so it can easily be positioned exactly where you want it for maximum scatter radiation protection.

- Raise and lower shield
- Smooth easy maneuvering
- Special spring mechanism allows for easy raising, lowering, and automatically holding shield
- .50mm LE double protective curtain
- 2 high-quality hospital grade locking casters



| Model Number | Description |
|----------------------|--|
| REV-MC | Economy Mobile Curtain |
| Dimensions | 66" H x 28" W |
| Curtain | 30" x 30" .50mm LE curtain with grey Revolution fabric |
| Curtain Color | Revolution Grey Custom colors available. |
| Frame Colors | White, Grey, or Black Custom colors available. Ships in standard white unless specified. |
| Frame Options | Powder Coated Steel Stainless Steel (Upgrade) |
| Wheels | 4 swivel casters, 2 with brakes |
| Crate Size | Box 1: 24" x 30" x 12" Box 2: 6" x 6" x 6" |
| Crate Weight | Box 1: 20 lbs / 9.07 kg Box 2: 27 lbs / 12.25 kg |

