


Mid-Tall Mobile Barrier

683473



- Powder coated steel frames, precision TIG welded for years of trouble free service completely enclosing the leaded glass panel for increased durability and safety
- Lead glass or lead acrylic options available
- 1.5 mm lead equivalent opaque base covered in "Sintra" plastic for easy cleaning and durability
- Modular parts can be easily replaced, including the glass if needed.
- Stainless steel finish available. Use -SS when ordering.
- Shipped in standard white unless specified.

Model Number	Description
683473	Mid-Tall Mobile Barrier - Lead
683473-AC	Mid-Tall Mobile Barrier - Acrylic
Dimensions	74½" H x 34½" W X 25" D
Glass Size	60" H x 30" W
Glass Options	Lead: 2.0 mm Pb LE Acrylic: .50 mm Pb LE
Frame Colors	White, Grey, or Black Custom colors available. Ships in standard white unless specified.
Frame Options	Powder Coated Steel Stainless Steel (Upgrade)
Curtain Color	Revolution Grey Custom colors available.
Wheels	4 swivel casters, 2 with brakes
Crate Size	73" x 30" x 43"
Crate Weight	380 lbs / 172.36 kg