


Physician Protection Mobile Screen

076994



A unique new design for operating physicians that allows unparalleled protection for close-to-table studies. Sturdy, yet thoughtfully engineered with a small footprint so it can easily be positioned exactly where you want it for maximum scatter radiation protection.

- Screen height adjustable by 11.8" (30 cm). Includes safety catch to set in desired position
- .50mm LE double protective curtain
- Horizontal position can be adjusted to fit close to table
- 2 high-quality hospital grade locking casters
- Powder coated steel frame
- Modular parts can be easily replaced, including the lead acrylic if needed.
- Shipping worldwide by ground, air or sea freight



Model Number	Description
076994	Physician Protection Mobile Screen
076994-IT	Optional 9"x 24" Stainless Steel Instrument Tray
Dimensions	Minimum Height: 57.5" (146 cm), Maximum Height: 69.7" (177 cm) Base: 27" W x 24" D
Glass Size	45.7" W x 18.7" H
Glass Options	.50mm LE acrylic
Curtain	Standard Revolution Grey Custom colors available.
Frame Colors	White, Grey, or Black Custom colors available. Ships in standard white unless specified.
Frame Options	Powder Coated Steel Stainless Steel (Upgrade)
Wheels	4 swivel casters, 2 with brakes
Crate Size	48" x 40" x 12"
Crate Weight	160 lbs / 72.57 kg

Optional 9"x 24" Stainless Steel Instrument Tray 076994-IT

