


Regular Mobile Barrier

683484



- Powder coated steel frames, precision TIG welded for years of trouble free service completely enclosing the leaded glass panel for increased durability and safety
- Lead glass or lead acrylic options available
- 1.5 mm lead equivalent opaque base covered in "Sintra" plastic for easy cleaning and durability
- Modular parts can be easily replaced, including the glass if needed.
- Stainless steel finish available. Use -SS when ordering.

| Model Number | Description |
|----------------------|--|
| 683484 | Regular Mobile Barrier - Lead |
| 683484-AC | Regular Mobile Barrier - Acrylic |
| Dimensions | 64½" H x 40½" W X 25" D |
| Glass Size | 24" H x 36" W |
| Glass Options | Lead: 2.0 mm Pb LE Acrylic: .50 mm Pb LE |
| Frame Colors | White, Grey, or Black Custom colors available. Ships in standard white unless specified. |
| Frame Options | Powder Coated Steel Stainless Steel (Upgrade) |
| Curtain Color | Revolution Grey Custom colors available. |
| Wheels | 4 swivel casters, 2 with brakes |
| Crate Size | 30" x 65" x 48" |
| Crate Weight | 430 lbs / 195.04 kg |