


Revolution Solid Panel Mobile Curtain

REV-SPMC



- Heavy-duty hospital-grade purpose-built unit.
- Adjustable height between 37" and 56"
- 30" wide by 30" tall lead vinyl curtain
- 2 high-quality hospital grade locking casters
- Smooth and easy maneuvering

Model Number	Description
REV-SPMC50	.50mm Pb Solid Panel Adjustable Height
REV-SPMC10	1.0mm Pb Solid Panel Adjustable Height
Dimensions	56" H x 30" W
Curtain Size	30" H X 30" W
Curtain Color	Revolution Grey Custom colors available.
Glass Options	Lead: 1.0 mm Pb LE Acrylic: .50 mm Pb LE
Frame Colors	White, Grey, or Black Custom colors available. Ships in standard white unless specified.
Frame Options	Powder Coated Steel Stainless Steel (Upgrade)
Wheels	4 swivel casters, 2 with brakes
Crate Size	Box 1: 21" x 21" x 10" Box 2: 16" x 33" x 7"
Crate Weight	Box 1: 18 lbs / 8.17 kg Box 2: 28 lbs / 12.7 kg



Adjustable height locking system from 37" to 56"