


“Shorty” Mobile Barrier

683458



- Powder coated steel frames, precision TIG welded for years of trouble free service completely enclosing the leaded glass panel for increased durability and safety
- Lead glass or lead acrylic options available
- 1.5 mm lead equivalent opaque base covered in “Sintra” plastic for easy cleaning and durability
- Modular parts can be easily replaced, including the glass if needed.
- Stainless steel finish available. Use -SS when ordering.

Model Number	Description
683458	Shorty Mobile Barrier - Lead
683458-AC	Shorty Mobile Barrier - Acrylic
Dimensions	45" H x 34½" W X 25" D
Glass Size	24" H x 30" W
Glass Options	Lead: 2.0 mm Pb LE Acrylic: .50 mm Pb LE
Frame Colors	White, Grey, or Black Custom colors available. Ships in standard white unless specified.
Frame Options	Powder Coated Steel Stainless Steel (Upgrade)
Curtain Color	Revolution Grey Custom colors available.
Wheels	4 swivel casters, 2 with brakes
Crate Size	45" x 30" x 42"
Crate Weight	220 lbs / 99.79 kg