

# Small Window Mobile Barrier

## 683465



- Powder coated steel frames, precision TIG welded for years of trouble free service completely enclosing the leaded glass panel for increased durability and safety
- Lead glass or lead acrylic options available
- 1.5 mm lead equivalent opaque base covered in "Sintra" plastic for easy cleaning and durability
- Modular parts can be easily replaced, including the glass if needed.
- Stainless steel finish available. Use -SS when ordering.
- Shipped in standard white unless specified

Model Number	Description
<b>683465</b>	Small Window Mobile Barrier - Lead
<b>683465-AC</b>	Small Window Mobile Barrier - Acrylic
<b>Dimensions</b>	72 1/4" H x 52 1/2" W X 25" D
<b>Glass Size</b>	12" H x 12" W
<b>Glass Options</b>	Lead: 2.0 mm Pb LE Acrylic: .50 mm Pb LE
<b>Frame Colors</b>	White, Grey, or Black Custom colors available. Ships in standard white unless specified.
<b>Frame Options</b>	Powder Coated Steel Stainless Steel (Upgrade)
<b>Floor Curtain Color</b>	Revolution Grey Custom colors available.
<b>Wheels</b>	4 swivel casters, 2 with brakes
<b>Crate Size</b>	30" x 72" x 60"
<b>Crate Weight</b>	380 lbs / 217.72 kg