


Technician Protection Mobile Screen

076993



The Technologist Protection X-ray Mobile Barrier can be easily raised and lowered by hand and stays in position thanks to the perfectly counter-balanced internal springs. The high quality steel construction and stainless springs provide a lifetime of trouble-free use.

- Retractable Window
- Multiple heights
- Perfectly counter balanced
- Easily raised and lowered
- 1.0mm Pb equivalent base
- .50mm Pb equivalent lead acrylic window
- 2 high-quality hospital grade locking casters
- Complete with bar to hold controls



Model Number	Description
076993	Technician Protection Mobile Screen
Size Raised	73" H x 31" W x 20 1/2" D
Size Lowered	44" H x 31" W x 20 1/2" D
Glass Size	37.4" H x 27" W
Glass Options	Lead: 2.0 mm Pb LE Acrylic: .50 mm Pb LE
Frame Colors	White, Grey, or Black Custom colors available. Ships in standard white unless specified.
Frame Options	Powder Coated Steel, Stainless Steel (Upgrade)
Wheels	4 swivel casters, 2 with brakes
Crate Size	48" x 40" x 10"
Crate Weight	180 lbs / 81.65 kg



New steel frame adds strength and prevents bending or warping of viewing window