


Wide Mobile Barrier

683482



- Powder coated steel frames, precision TIG welded for years of trouble free service completely enclosing the leaded glass panel for increased durability and safety
- Lead glass or lead acrylic options available
- 1.5 mm lead equivalent opaque base covered in "Sintra" plastic for easy cleaning and durability
- Modular parts can be easily replaced, including the glass if needed.
- Stainless steel finish available. Use -SS when ordering.
- Shipped in standard white unless specified.

Model Number	Description
683482	Wide Mobile Barrier - Lead
683482-AC	Wide Mobile Barrier - Acrylic
Dimensions	72" H x 52½" W X 25" D
Weight	235 lbs.
Window	36" H x 48" W
Glass	Lead or Acrylic (-AC)
Frame Colors	White, Grey, or Black Custom colors available. Ships in standard white unless specified.
Frame Options	Powder Coated Steel Stainless Steel (Upgrade)
Curtain Color	Revolution Grey Custom colors available.
Wheels	4 swivel casters, 2 with brakes
Crate Size	30" x 72" x 60"
Crate Weight	480 lbs / 217.72 kg