

# Achilles Lead Glasses



**800076**  
Blue



**800075**  
Red



**800077**  
Plated Silver



**800079**  
Shiny Black

Lens Width: 53mm



Bridge Width: 17mm



Lens Height: 35mm



Frame Width: 128mm



- Vented radiation protection goggle with neoprene foam on the inside designed to protect against fluid splashing and perspiration
- Easy to adjust with elastic head band
- Prescription compatible (for most Rx)
- Lightweight and comfortable
- Large lens area for maximum protection

<b>Frame Size</b>	Larger
<b>Prescription</b>	Plano Single Vision (-SV) Bifocal (-BF) Bifocal Progressive (-BF-P)
<b>Colors</b>	800076 Blue 800075 Red 800077 Plated Silver 800079 Shiny Black
<b>Weight</b>	69G
<b>Protection</b>	.75 mm LE
<b>Includes</b>	Anti-fog spray, safety strap add-on, protective case