

Cyberflex Lead Glasses Single Vision, Plano



683206-S Clear/Blue

Lens Width: 59mm

Bridge Width: 19mm

Frame Width: 136mm



Temple: 135-150mm

Lens Height: 40mm



- High Quality Safety Glasses
- Adjustable temple bar
- Large viewing area
- Side Protection, 0.50mm
- Oversized protective lens area

Frame Size	Universal
Prescription	Single Vision (-SV) Plano
Colors	Clear/Blue
Weight	95.4 G
Protection	.75 mm LE
Includes	Anti-fog spray, safety strap add-on, protective case