

Cyberflex Lead Glasses



683206 Clear/Blue

Lens Width: 59mm

Bridge Width: 19mm

Frame Width: 136mm



Temple: 135-150mm

Lens Height: 40mm



- Optional .50mm LE side shields
- Adjustable temple bar
- Large viewing area
- Built-in side shields for added splash protection
- Oversized protective lens area
- Splash protection

Frame Size	Universal
Prescription	Single Vision (-SV) Plano
Colors	683206 Clear/Blue
Weight	83G
Protection	.75 mm LE
Includes	Anti-fog spray, safety strap add-on, protective case