

Fitovers EcoFit



832001
Tortoise



832000
Black

Fits over Rx glasses
not exceeding width:
135mm



Fits over Rx glasses
not exceeding height:
41mm



- Fits frames not exceeding 153mm width x 41mm height
- Wrap around Radiation protection
- .50 mm side shield protection
- Hypo-allergenic Crystal Nylon™ Frame
- Designed to be worn over prescription eyewear

Frame Size	Universal
Prescription	Plano
Colors	Black, Tortoise
Weight	89G
Protection	.75 mm LE
Side Protection	.50 mm LE
Includes	Side shields, Anti-fog spray, safety strap add-on, protective case