

## Fitovers Lotus



**800304**  
Matte Black



**800306**  
Smoke Marble



**800307**  
Claret Stripe

- Fits frames not exceeding 141mm width x 40mm height.
- Wrap around Radiation protection
- .50 mm side shield protection
- Hypo-allergenic Crystal Nylon™ Frame
- Designed to be worn over medium, oval, or round prescription eyewear frame shapes

Fits over Rx glasses  
not exceeding width:  
**141mm**



Fits over Rx glasses  
not exceeding height:  
**40mm**



Frame Size	Universal
Prescription	Plano
Colors	Matte Black, Smoke Marble, Claret Stripe
Weight	83G
Protection	.75 mm LE
Side Protection	.50 mm LE
Includes	Anti-fog spray, safety strap add-on, protective case