

Fitovers Orion



ORION-01
Black



ORION-02
Blue Marble

- Fits frames not exceeding 140mm width x 43mm height
- .50 mm side shield protection
- Wrap around radiation protection
- Hypo-allergenic Crystal Nylon™ Frame
- Designed to be worn over medium prescription eyewear frame shapes
- Non Rx

Fits over Rx glasses
not exceeding width:
140mm



Fits over Rx glasses
not exceeding height:
43mm



Frame Size	Universal
Prescription	Plano
Colors	Black, Blue Marble
Weight	84G
Protection	.75 mm LE
Side Protection	.50 mm LE
Includes	Anti-fog spray, safety strap add-on, protective case