

GUCCI

RADIATION PROTECTION GLASSES



GG0560-005
Black / Crystal


Lens Width: 55mm




Bridge Width: 20mm



Temple Length: 145mm



Lens Height: 41.6mm



Anima Décor inspired by street styles trend for a young, cool and fresh look. Masculine squared acetate frame embellished with transparent temples which highlight the internal wire core with the Gucci graphic band logo. Classic black tones alternated with transparent crystal ones.

Frame Size	55mm - 20mm - 145mm
Prescription	Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive
Colors	Black/Crystal
Weight	73 G
Protection	.75 mm LE
Includes	Case, Cloth & Anti-Fog Spray