

GUCCI

RADIATION PROTECTION GLASSES



GG0758-003
Grey


Lens Width: 56mm




Bridge Width: 17mm



Temple Length: 150mm



Lens Height: 38.6mm



Best Seller shape interpreted in a modern, light and high quality version. Rectangular masculine shape. Temples with Gucci GG lettering logo. Classic color palette with solid and transparent tones.

Frame Size	56mm - 17mm - 150mm Alternative Fit
Prescription	Plano, Single Vision, Most Bi-Focal & Bi-Focal Progressive
Colors	Grey
Weight	65 G
Protection	.75 mm LE
Includes	Case, Cloth & Anti-Fog Spray