

Nike 7082 Lead Glasses

Plano



7082-01
Matte Midnight Navy - Racer Blue



7082-02
Matte Gridiron - Gym Red



7082-03
Matte Sequoia - Volt

Lens Width: 57mm



Bridge Width: 16mm



Temple Length: 145mm



Lens Height: 37.9mm



- Universal fit nosepiece with two replaceable nosebridge sizings for universal fitting.
- Rubber nosepiece and temple tips designed to become tackier as the skin perspires for a strong and secure eyewear fit.
- Core wire memory metal temple tips for on the fly adjustability and grip.
- Lightweight, durable nylon frame

FRAME SIZE:	Universal
WEIGHT:	64.4G
PRESCRIPTION:	Plano
COLORS:	Matte Midnight Navy - Racer Blue, Matte Gridiron - Gym Red, Matte Sequoia - Volt
PROTECTION	.75mm LE
INCLUDES	Anti-fog spray, safety strap add-on, protective case