



Nike 7084 Lead Glasses



7084-01
Matte Black



7084-02
Matte Anthracite



7084-03
Matte Space Blue

Bridge Width: 15mm



Frame Width: 127mm

Temple Length: 140mm



Lens Width: 53.9mm



Lens Height: 35.8mm



- Injected TR90, lightweight frame front
- Dual injected nylon end piece and cam hinges
- Rubber temple sleeves with core wire for adjustability
- Dual injected rubber and nylon, snap on ventilated nose pad to reduce fogging
- Alternative fit nose bridge included for wider range of fit
- Full-Rim

Frame Size	Universal
Prescription	Single Vision (-SV) Plano
Colors	7084-01 Matte Black 7084-02 Matte Anthracite 7084-03 Matte Space Blue
Weight	55.4 G
Protection	.75 mm LE
Includes	Anti-fog spray, safety strap add-on, protective case