



# Nike 7101 Lead Glasses



**7101-01**  
Matte Anthracite

Lens Width: 53mm



Bridge Width: 15mm



Temple Width: 140mm



Lens Width: 37.1mm



Nike 7101 is a full rim frame made of plastic. It features the Nike swoosh logo on temples. This product offers our standard high quality, distortion-free radiation reduction lenses.

<b>Prescription</b>	Single Vision (-SV) Plano
<b>Frame Size</b>	Universal
<b>Colors</b>	7101-01 Matte Anthracite
<b>Weight</b>	58G
<b>Protection</b>	.75 mm LE
<b>Includes</b>	Anti-fog spray, safety strap add-on, protective case