



# Nike 7237 Lead Glasses



**7237-01**  
Matte Black-Dark Grey

Bridge: 17mm



Arm: 140mm



Lens: 52mm W



Lens: 31 mm H



- Lightweight acetate frame design for comfort and durability
- Stacked acetate temples with pop color accents in sneaker tread pattern
- Cam action hinges: snap firmly into place
- Jeweled plaque Nike Swoosh logo for pop color accent

<b>Frame Size</b>	Universal
<b>Prescription</b>	Single Vision (-SV) Plano
<b>Colors</b>	7237 - Matte Black-Dark Grey
<b>Weight</b>	56G
<b>Protection</b>	.75 mm LE
<b>Includes</b>	Anti-fog spray, safety strap add-on, protective case