

## Nike 7116 Lead Glasses

Single Vision, Plano



**7116-01**  
Matte Black - Volt



**7116-02**  
Matte Dark Grey - Black



**7116-03**  
Matte Gridiron - Pacific Blue

**Bridge Width:** 15mm



**Temple:** 140mm



**Lens Width:** 56.1mm



**Lens Height:** 38.8mm



- FRAME SIZE:** Universal
- WEIGHT:** 65.1G
- PRESCRIPTION:** Single Vision, Plano
- COLORS:** Matte Black/Grey, Matte Gery/Blue
- PROTECTION:** .75mm LE
- INCLUDES:** Anti-fog spray, safety strap add-on, protective case

- Rubber nosepad and temple tip inlays designed to become tackier as the skin perspires for a strong and secure eyewear fit.
- Core wire memory metal temple tips for on the fly eyewear adjustability and grip.
- Lightweight, durable nylon frame