

Nike 7117 Lead Glasses

Plano



7117-01
Matte Black - Gym Red



7117-02
Matte Midnight Navy - Obsidian

Lens Width: 54mm



Lens Height: 38.9mm



Bridge Width: 16mm



Temple Length: 140mm



- Rubber nosepad and temple tip inlays designed to become tackier as the skin perspires for a strong and secure eyewear fit.
- Core wire memory metal temple tips for on the fly eyewear adjustability and grip.
- Lightweight, durable nylon frame

FRAME SIZE: Universal

WEIGHT: 60G

PRESCRIPTION: Plano

COLORS: Matte Black - Gym Red, Matte Midnight Navy - Obsidian

PROTECTION: .75mm LE

INCLUDES: Anti-fog spray, safety strap add-on, protective case