

Nike 7238 Lead Glasses

Single Vision, Plano



7238-04
Matte Black/Grey



7238-05
Matte Grey/Blue

Bridge Width: 18mm



Temple: 140mm



Lens Width: 52mm



Lens Height: 36.4mm



- FRAME SIZE:** Universal
- WEIGHT:** 58 G
- PRESCRIPTION:** Single Vision, Plano
- COLORS:** Matte Black/Grey, Matte Grey/Blue
- PROTECTION:** .75mm LE
- INCLUDES:** Anti-fog spray, safety strap add-on, protective case

- Lightweight acetate frame design for comfort and durability.
- Stacked acetate temples with pop color accents in tread pattern
- Cam action hinges snap firmly into place
- Jeweled plaque Nike Swoosh logo for pop color accent