



Nike Bandit Lead Glasses



BANDIT-01
Black



BANDIT-02
Blue

The Nike Bandit is a full rim square frame which is made of plastic. It comes with nose pads that snap into the pad box. This product offers our standard high quality, distortion-free radiation-reducing lenses with .75mm lead equivalency.

- Ventilated nose bridge for reduced fogging and secure grip
- Tiny holes in top of frame for air circulation and fog/ moisture prevention



Lens Width: 59mm



Bridge Width: 17mm



Temple Length: 145mm



Lens Height: 45mm



Frame Size	Large Frame
Prescription	Plano
Colors	BANDIT-01 Black BANDIT-02 Blue
Weight	80G
Protection	.75 mm LE
Includes	Anti-fog spray, safety strap add-on, protective case