



Nike Brazen Fury Lead Glasses

Plano



BRAZENFURY-01
Matte Dark Grey / Silver Flash



BRAZENFURY-02
Matte Black / White Red



BRAZENFURY-03
Matte Midnight Navy

Lens Width: 55mm



Lens Height: 36.5mm



Bridge Width: 14mm



Temple Length: 130mm



- Rubber nose pad and temple tip inlays designed to become tackier as the skin perspires for a strong and secure eyewear fit.
- Core wire memory metal temple tips for on the fly eyewear adjustability and grip.
- Lightweight, durable nylon frame

FRAME SIZE: Universal

PRESCRIPTION: Plano

COLORS: Matte Dark Grey / Silver Flash, Matte Black / White Red, Matte Midnight Navy

PROTECTION: .75mm LE

INCLUDES: Anti-fog spray, safety strap add-on, protective case