



Nike Brazen Lead Glasses



BRAZEN-02
Matte Black / Military Blue



BRAZEN-04
Shiny Black

Lens Width: 60mm



Bridge Width: 16mm



Frame Width: 135mm



Lens Height: 40mm



The Nike Brazen glasses combine streamlined style, comfort and lightweight durability. They feature a nylon frame and ventilated nose bridge. The Brazen is our most popular model.

Frame Size	Universal
Prescription	Plano
Colors	BRAZEN-02 Matte Black / Military Blue BRAZEN-03 Shiny Black
Weight	77G
Protection	.75 mm LE
Includes	Anti-fog spray, safety strap add-on, protective case