




Nike Current Lead Glasses



CURRENT-01
Matte Black

Bridge Width: 19mm



Frame Width: 127mm

Temple Length: 145mm



Lens Width: 52mm



Lens Height: 45.7mm



Nike Current sunglasses combine everyday appeal with performance benefits. With a nod to the aviator style, they provide clear vision and a stable fit in a super-low profile frame.

Frame Size	Universal
Prescription	Single Vision (-SV) Plano
Colors	Current-01 Matte Black
Weight	60.2G
Protection	.75 mm LE
Includes	Anti-fog spray, safety strap add-on, protective case