



# Nike Fleet Lead Glasses



**FLEET-01**  
Matte Black/Volt



**FLEET-03**  
Matte Blue/Grey

Lens Width: 55mm



Bridge Width: 16mm



Temple Length: 135mm



Lens Height: 42mm



The Nike Fleet is a full rim square frame which is made of plastic. It comes with nose pads that snap into the pad box. This product offers our standard high quality, distortion-free radiation-reducing lenses with .75mm lead equivalency. This model is a smaller version of the Nike Bandit.

<b>Frame Size</b>	Universal
<b>Prescription</b>	Plano
<b>Colors</b>	FLEET-01 Matte Black/Volt FLEET-03 Matte Blue/Grey
<b>Weight</b>	71G
<b>Protection</b>	.75 mm LE
<b>Includes</b>	Anti-fog spray, safety strap add-on, protective case