



Nike Ledge Lead Glasses



LEDGE-01
Anthracite

Bridge Width: 19mm



Frame Width: 144mm

Arm Length: 145mm



Lens Width: 56mm



Lens Height: 42.8mm



With a super-low profile and strong, squared off design, the Nike Ledge lead glasses provide style and comfort.

Frame Size	Universal
Prescription	Single Vision (-SV) Plano
Colors	LEDGE-01 Matte Black
Weight	65.1G
Protection	.75 mm LE
Includes	Anti-fog spray, safety strap add-on, protective case