



Nike Maverick RGE Lead Glasses

Plano



MAVERICKRGE-01
Matte Black / Green



MAVERICKRGE-02
Matte Black / Grey



MAVERICKRGE-03
Matte Midnight Navy / Silver Flash

Lens Width: 59mm



Lens Height: 40.4mm



Bridge Width: 15mm



Temple Length: 145mm



Full rubber nosebridge and temple tips designed to become tackier as the skin perspires for a strong and secure eyewear fit.

- FRAME SIZE:** Universal
- PRESCRIPTION:** Plano
- COLORS:** Matte Black / Green
- PROTECTION:** .75mm LE
- INCLUDES:** Anti-fog spray, safety strap add-on, protective case