



Nike Premier Lead Glasses



PREMIER-04
Cool Grey



PREMIER-05
Matte Blue

Bridge Width: 13mm



Frame Width: 145mm



Lens Width: 60mm



Temple Length: 135mm



Lens Height: 39.3mm

Lightweight and comfortable, the Nike Premier lead glasses feature slimmed down temples with tips that grip for stability.

Frame Size	Universal
Prescription	Plano
Colors	PREMIER-04 Nike Premier Cool Grey PREMIER-05 Nike Premier Matte Blue
Dimensions	60-13-135 B:39.3 ED:65.2 FW:145
Weight	71.8G
Protection	.75 mm LE
Includes	Anti-fog spray, safety strap add-on, protective case